

The field of drug testing is constantly changing. With many years of experience in developing assays, we are able to anticipate future drug testing needs and to provide innovative total test solutions. Offering the most extensive menu in the industry we continue helping our customers to meet their distinct drug testing requirements.

Thermo Fisher Scientific


CDx 90



### The SMALL & SMART dedicated drug testing solution

The CDx 90 Analyser is especially developed to meet the requirements of low volume, dedicated drug testing laboratories.

The CDx 90 is SMALL, so it can flawlessly fit the low volume testing lab needs in terms of space and throughput. It is as SMART as "big" analysers and has all necessary features and trappings that ensure accurate results, high-level convenience and ease-of-use. A complete drug testing solution from one supplier, a real "one-stop shop".



### Key Features:

#### Drug Testing Solution

- Compact and convenient benchtop analyser
- Innovative, unique and full menu of Microgenics reagents
- DAT, Validity, Tox, TDM and ISD testing – reagents, calibrators & controls
- Reliable and consistent results for the sample of your choice
- Intuitive system software
- Flexible data management tool (Toxicom® optional)
- Comprehensive instrument and product training
- Responsible scientific support and service

#### Analyser

- Fully automated, random access analyser
- Throughput – 90 tests / hour
- STAT sample – interruption permitted, even during analysis
- Up to 20 reagent containers and up to 10 samples on-board
- Reagent and sample cooling
- Reagent and sample liquid level sensing
- Reagent and sample external positive barcode identification
- Touchscreen
- Low Maintenance
- 10 programmable open channels

## Thermo Fisher CDx 90 – SMALL & SMART

The analyser runs the comprehensive range of Microgenics Drugs of Abuse, Toxicology, Therapeutic Drug Monitoring and Immunosuppressive Drug Testing assays (for the full menu see next page), including several unique parameters, e.g.:



### Ethyl Glucuronide (EtG)

The first fully automated immunoassay for the detection of EtG in urine. EtG has proven to be an excellent biomarker for recent alcohol consumption; it has a significantly longer detection time than Alcohol (80 vs 8 hrs).

### High Sensitive Benzodiazepine

This CEDIA<sup>®</sup> assay combines:

- Integrated online hydrolysis for the detection of glucuronidated Benzodiazepines (e.g. Oxazepam and Temazepam)
- The recognition of 7-amino-metabolites from the new clinical important and abused Benzodiazepines (e.g. Flunitrazepam and Clonazepam)

### Methadone Metabolite (EDDP)

Testing for the Methadone Metabolite EDDP provides conclusive evidence of Methadone ingestion:

- Identification of specimens that have been spiked with Methadone;
- Renal clearance of EDDP is unaffected by urinary pH resulting in fewer false negatives;
- Protection of patients with rapid Methadone metabolism or on low dose regimens

### Lithium


The only liquid ready to use and fully automated Lithium assay available to guarantee fast turnaround of excellent quality results.

### Antiepileptic Drug Assays

Besides the common Antiepileptic Drug Assays like Carbamazepine and Phenytoin, Microgenics has also developed excellent assays for the newer Antiepileptic drugs like Lamotrigine and Topiramate.


### MAS<sup>®</sup> Quality Control Programme

This QC Solution can be used in combination with the CDx 90 and the Microgenics drug assays. The programme consists of Control Material and LabLink<sup>®</sup> xL, a web-based quality assurance programme.





## Instrumentation



### Thermo Fisher CDx 90 – SMALL & SMART – Test Menu

#### Drugs of Abuse Assays

Amphetamines / Ecstasy  
 Barbiturates  
 Benzodiazepines  
 Buprenorphine  
 Cannabinoids (multi level)  
 Cocaine Metabolite  
 Cotinine  
 Ecstasy  
 Ethyl Alcohol  
 Ethyl Glucuronide  
 Heroin Metabolite (6-AM)  
 Ketamine  
 LSD  
 Methadone  
 Methadone Metabolite (EDDP)  
 Methaqualone  
 Opiates  
 Oxycodone  
 Phencyclidine  
 Propoxyphene

#### Sample Validity Assays

CEDIA® Sample Check  
 Creatinine

#### Alternative Samples

Drugs of Abuse - Oral Fluid  
 Drugs of Abuse - Hair

#### Serum Toxicology Assays

Acetaminophen  
 Barbiturates  
 Benzodiazepines  
 Salicylate  
 Tricyclic Antidepressants

#### Therapeutic Drug Monitoring

##### *Antiepileptic Drug Assays*

Carbamazepine  
 Lamotrigine  
 Phenobarbital  
 Phenytoin  
 Topiramate  
 Valproic Acid  
 Zonisamide

##### *Antimicrobial Drug Assays*

Amikacin  
 Gentamicin  
 Teicoplanin  
 Tobramycin  
 Vancomycin

##### *Immunosuppressive Drug Assays*

Cyclosporine (Low & High)  
 MPA

##### *Cardioactive Drug Assays*

Digitoxin  
 Digoxin  
 Procainamide  
 N-Acetylprocainamide (NAPA)  
 Quinidine

##### *Antiasthmatic Drug Assay*

Theophylline

##### *Antidepressant Drug Assay*

Lithium

## Thermo Fisher CDx 90 - SMALL & SMART Analyser Specifications

### Technical Data

#### General

<b>System</b>	Fully automated, random access benchtop analyser
<b>Throughput</b>	90 tests per hour
<b>Category</b>	Dedicated drug testing system (photometric)
<b>Assay menu</b>	Drugs of Abuse Testing (DAT - urine) Validity testing (Adulteration - urine) Toxicology testing (Tox - serum) Therapeutic Drug Monitoring (TDM) Immunosuppressive Drug Monitoring (ISD)

#### Sample Handling

<b>Sample type</b>	Whole blood, serum, plasma, urine, oral fluid (saliva), hair and others
<b>Sample container</b>	Cups and primary tubes (diameter 13 - 16 mm, length 75 - 100 mm)
<b>Sample input</b>	Removable tray with sample cup/tube holder on turntable
<b>Sample cooling</b>	Peltier element cooling (8 - 15° Celsius)
<b>Sample capacity</b>	Up to 10 samples, cups and/or tubes
<b>STAT sample</b>	Interruption permitted, even during analysis
<b>Sample volume</b>	2 to 35 µL (0.1 µL/step)
<b>Sample dilution</b>	Automatic, programmable sample dilution
<b>Sample identification</b>	Position number and/or Positive sample identification with handheld barcode reader
<b>Sample probe</b>	Micropipette / microsyringe (2 to 35 µL - 0.1 µL/step) Liquid level sensor, capacity measurement Washed inside and outside with purified water and/or detergent(s)

#### Reagent Handling

<b>Reagent type</b>	CEDIA®, DRI®, QMS® and Infinity™ drug assays
<b>Reagent container</b>	20 and 50 ml reagent bottles
<b>Reagent input</b>	Removable tray with reagent container on turntable, accommodates Microgenics drug testing reagent bottles
<b>Reagent cooling</b>	Peltier element cooling (8 - 15° Celsius)
<b>Reagent capacity</b>	Up to 20 reagent bottles: 10 x 20 mL / 10 x 50 mL, or with adapters 20 x 20 mL
<b>Reagent volume</b>	20 to 350 µL (1.0 µL/step)
<b>Inventory</b>	Calculation of remaining reagent volume
<b>Reagent identification</b>	Position number and/or Positive reagent identification with handheld barcode reader
<b>Reagent probe</b>	Micropipette / microsyringe (20 to 350 µL - 1.0 µL/step) Liquid level sensor, capacity measurement Washed inside and outside with purified water and/or detergent(s)

#### Reaction

<b>Cuvettes</b>	Disposable resin cuvettes
<b>Cuvette ring</b>	Maximum 24 cuvettes
<b>Cuvette pack</b>	96 single cuvettes per tray
<b>Cuvette volume</b>	Minimum 120 µL, maximum 450 µL
<b>Cuvette temperature</b>	Direct heating, 37°C (+/- 0.3°C)
<b>Cuvette handling</b>	Automatic feeding and disposal
<b>Pipetting</b>	1 pipettor for sample and reagents
<b>Dispensing</b>	Sample 100 µL syringe, reagent 1.000 µL syringe
<b>Mixing method</b>	Stirring paddle
<b>Incubation time</b>	Five minutes after adding sample for one reagent assay Five minutes after adding second reagent for two reagent assay

#### Measurement

<b>Method</b>	Photometric, absorbance measurement Bi-chromatic and mono-chromatic
<b>Wavelength</b>	8 wavelengths (filter wheel 340 – 800 nm)
<b>Light source</b>	Halogen tungsten lamp
<b>Linear range</b>	OD 0.0 – 3.5
<b>Analytical method</b>	Endpoint and Rate

#### User Interface

<b>Interface</b>	Built-in PC with touchscreen and handheld barcode reader
<b>USB connection</b>	Mouse, keyboard, printer and memory stick(s)
<b>Status display</b>	Graphic and colour coded operation status display
<b>Calibration monitor</b>	Factor, linear, logit/log, spline & exponential Numeric and graphic display
<b>Quality control</b>	Westgard multi-rule - numeric and graphic display
<b>Reaction monitor</b>	Reaction curve - numeric and graphic display

#### Installation

<b>Dimension</b>	620 mm (W) x 555 mm (D) x 570 mm (H)
<b>Weight</b>	65 kg
<b>Optional</b>	Handheld barcode reader Water and waste liquid sensors Data management tool Toxicom®